

JohnsonDiversey ClearKlens Tego 2000 SC

Terminal disinfectant for open plant applications

VH25

Description

ClearKlens Tego® 2000 SC is a concentrated disinfectant for the Pharmaceutical industry.

Key properties

- **ClearKlens Tego® 2000 SC** is a leading terminal disinfectant in the Pharmaceutical industry.
- **ClearKlens Tego® 2000 SC** is both 0.2µm filtered and gamma irradiated.
- **ClearKlens Tego® 2000 SC** is triple bagged for hygiene security.
- **ClearKlens Tego® 2000 SC** can be applied as a spray or by means of a clean cloth. Use for the disinfection of walls, floors, work surfaces and process equipment.

Small items or dismantled equipment may be immersed in a **ClearKlens Tego® 2000 SC** solution.

- **ClearKlens Tego® 2000 SC** is based on an amphoteric surfactant system.
- **ClearKlens Tego® 2000 SC** is a 50ml single shot triple bagged disinfectant, which has been filtered through a 0.2µm filter and irradiated.
- **ClearKlens Tego® 2000 SC** can also be used for fogging.

ClearKlens Tego® 2000 is the registered trademark of Degussa AG, Germany.

Benefits

- Sterile pack
- Detex labelled for proof of irradiation
- Triple bagged for easy transfer into the cleanroom
- Odourless
- Non tainting
- Low toxicity
- Broad spectrum of activity
- Effective in soft or hard water

Use instructions

ClearKlens Tego® 2000 SC is usually used at 1% v/v aqueous solution (empty the contents of a 50ml bottle into 5 litres of water)

For fogging use **ClearKlens Tego® 2000 SC** as a 2% solution at approximately 1L per 100m³.

Technical data

Appearance Clear:	Colourless to pale yellow liquid
Relative density at 20°C	0.998
pH (neat solution at 20°C)	8.0
Chemical Oxygen Demand (COD)	425mg O ₂ /g
Biochemical Oxygen Demand (BOD)	0.142mg O ₂ /g
Nitrogen Content (N)	19 g/Kg
Phosphorous Content (P)	none

The above data is typical of normal production and should not be taken as a specification.



JohnsonDiversey ClearKlens Tego 2000 SC

Safe handling and storage information

Store in original closed containers, away from sunlight and extremes of temperature. Full guidance on the handling and disposal of this product is provided in a separate Material Safety Data Sheet.
For professional use only

Product compatibility

ClearKlens Tego® 2000 SC when applied at the recommended concentration and temperature is suitable for use on the grades of stainless steel commonly found in the Pharmaceutical industry. In the event of uncertainty it is advisable to evaluate individual materials before any prolonged use.

Test method

Validated Titration Test Method available

Microbiological data

ClearKlens Tego® 2000 SC has been extensively tested using CEN 1276 Bactericidal Suspension Test. This has shown **ClearKlens Tego® 2000 SC** to have a broad spectrum of activity against both Gram-positive and Gram-negative bacteria.
Efficacy data available on request.

Available pack sizes

Article code	Pack size	Pack type
7513387	50 × 50ml	

JohnsonDiversey (UK)
Weston Favell Centre
Northampton NN3 8PD
Tel. 0800 525 525
www.johnsondiversey.com