

JohnsonDiversey ClearKlens Tego 2001 SS

ClearKlens Tego® 2001SS VH24

VH24

Description

ClearKlens Tego® 2001 SS is a ready to use disinfectant for the cleanroom industry. **ClearKlens Tego® 2001 SS** is a surface active, antimicrobial Amphoteric system, based on 1-alkyl-1,5-diazapentane in aqueous solution. **ClearKlens Tego® 2001 SS** is a 900ml bag in bottle ready to use solution, double bagged disinfectant, which has been filtered to 0.2µm and irradiated.

Key properties

- **ClearKlens Tego® 2001 SS** is the leading terminal disinfectant in the Pharmaceutical cleanroom industry.
 - **ClearKlens Tego® 2001 SS** is a ready to use solution in a bag in bottle delivery system to eliminated product contamination during use.
 - **ClearKlens Tego® 2001 SS** is both 0.2µm filtered and gamma irradiated.
 - ClearKlens Tego® 2001 SS is double bagged.
 - **ClearKlens Tego® 2001 SS** is used for the disinfection of work surfaces and process equipment.
 - **ClearKlens Tego® 2001 SS** should be applied as a spray.
- Tego® is a registered trademark of Degussa

Benefits

- 0.2µm filtered
- Gamma Irradiated
- Detex labelled for proof of irradiation
- Double bagged for easy transfer into the cleanroom
- Odourless
- Non tainting
- Low toxicity
- Broad spectrum of activity

Use instructions

ClearKlens Tego® 2001 SS is designed for spraying and wiping applications. Hold spray approximately 30 - 40 cm from area to be treated. Spray directly onto surface to ensure complete coverage. Allow sufficient wet contact time before proceeding. Always close the nozzle after use.

Technical data

Appearance:	Pale yellow liquid
Relative density at 20°C:	1.002
pH (neat solution at 20°C):	7.5
Chemical Oxygen Demand (COD)	9 mg O ₂ /g

The above data is typical of normal production and should not be taken as a specification.



JohnsonDiversey ClearKlens Tego 2001 SS

Safe handling and storage information

Store in original closed containers away from sunlight and extremes of temperature.
Full guidance on the handling and disposal of this product is provided in a separate Material Safety Data Sheet.
For professional users only

Product compatibility

Applications solutions, when used as directed, will not affect materials normally encountered in the aforementioned industries.

Test method

Validated Titration Test Method

Microbiological data

ClearKlens Tego 2001 SS has been extensively tested using CEN 1276 Bactericidal Suspension Test, the CEN 1650 Fungicidal Activity Suspension Test and the CEN 13697 Surface Test the new industry standard method. These have shown **ClearKlens Tego 2001 SS** to have a broad spectrum of activity against both Gram-positive and Gram-negative bacteria, yeasts and fungi.
Efficacy data is available on request

Available pack sizes

Article code	Pack size
7513403	6 × 0.9L

JohnsonDiversey (UK)
Weston Favell Centre
Northampton NN3 8PD
Tel. 0800 525 525
www.johnsondiversevy.com